

# Five Scanning Tips

---

- 1.** Use the **waist belt** provided while scanning and when taking length measurements: This helps ensure a consistent reference point.
- 2.** Scan patient wearing a **stockinette**: This will result in a higher quality scan.
- 3.** Take **measurements**: All appropriate measurements for us to compare to the scan.
- 4.** Take a standing, **one-piece scan** with knees bent slightly to reduce lordosis.
- 5.** Take **photos**: Use the SpinalTech3D™ Scan App and attach a couple photos of the patient along with submitting the scan.

