

Radiology Protocol



Take a Successful Providence Nocturnal Scoliosis® Orthosis In-Brace X-Ray

In-Brace X-Rays are critical for evaluating the effectiveness of any scoliosis orthosis. These X-Rays are the primary tool for determining the impact the brace is having on the patient internally and they help to inform the Orthotist of any adjustments or modifications that might be needed for improving correction and outcomes.

Being able to review In-Brace X-Rays, regardless of the brace's style, is important in evaluating mold modifications and ensuring effective brace blueprint design.

This will also help our team of scoliosis specialists successfully maintain the highest quality for you (the clinician), your scoliosis referrals, and, most importantly, your patient.

Providence Nocturnal Scoliosis® Orthosis



Out of Brace



In-Brace (supine)



IMPORTANT



To properly evaluate the design and effectiveness of the orthosis:

- The patient must be wearing the brace positioned on their body exactly as they have been instructed to do so by their Orthotist.
- The brace must be worn at full tightness.
- **Providence® Brace X-Rays must be taken with the patient lying down in supine position**

Email your In-Brace X-Rays to: inbracexrays@spinaltech.com

SpinalTechnology.com

Follow us:



**Spinal
Technology**



📍 **191 Mid-Tech Drive
West Yarmouth, MA 02673**

☎ **800.253.7868**