



Anterior View

Product Overview

The Kyphologic style w/ Anterior opening & "Cow Horns" that attach laterally and allows for a corrective force at the apex of kyphosis. The brace is used to correct Kyphosis and stabilize adult patients.

Spinal Indications

- Mild Structural Kyphosis
- Postural Stabilization
- Scheuermann's Kyphosis

Features

- The trims are designed to allow hips & scapulae to be relieved of pressure, to move more freely and is easier to conceal under clothing.
- "Cow Horns" allow for adjustment of height, & pivot for more extension force.



Posterior View



Left Lateral View



Right Lateral View

Send your orders to:
orders@spinaltech.com

Send your CAD files to:
cad@spinaltech.com

**For product options see
reverse side.**

Product Accessories



Posterior Left Lateral View

Posterior Right Lateral View



Anterior View with Channels

Recommended Plastic

Copolymer (CP) 1/8" or 5/32" depending on patient size
(Call for available color options)

Available Plastics

Polyethylene (PE) 3/32", 1/8", 5/32", 3/16" (Least Rigid)
Modified Polyethylene (MPE) 1/8", 5/32", 3/16" (Semi Rigid)

Recommended Aliplast Liner

3/16" Aliplast unlined

Other Liner Thickness Available

1/8", 1/4" Aliplast
Note: This orthosis can also be unlined.

Sternal attachments

Cow Horns
Cash Pad

Fasteners

Screws and Barrel Nuts
Speedy Rivets

ST 3.26