

Custom Spinal LSO

Anterior Overlap Orthosis

With Tongue



Anterior View



Left Lateral
View



Posterior View

Product Overview

This LSO utilizes a single opening. One piece of plastic is measured to wrap around the patient and tighten in the front (anterior), using three straps. An additional piece of plastic (tongue) is placed behind the straps. This brace is a good solution for patients that need total circumferential control.

Spinal Indications

- Post-Operative Lumbar Surgery
- Degenerative Disc Disease
- Spondylolysis
- Spinal Stenosis
- Compression Fracture
- Herniated Disc
- Low Back Pain

Features

- More total circumferential control.
- More total circumferential control.

Opening Available

Anterior Opening with Tongue (Pictured), Posterior Opening, Lateral Opening (Right or Left)

Send your orders to:
orders@spinaltech.com

Send your CAD files to:
cad@spinaltech.com

For product options see
reverse side.

Product Options

Recommended Plastic

Polyethylene (PE) 3/32", 1/8", 5/32", 3/16" (Least Rigid)

Modified Polyethylene (MPE) 1/8", 5/32", 3/16" (Semi Rigid) (if Transfer pattern is requested)

Available Plastics

Copolymer (CP) 1/8" or 5/32" depending on patient size
(Call for available color options)

Recommended Aliplast Liner

3/16" Aliplast

Other Liner Thickness Available

1/8", 1/4" Aliplast

Note: This orthosis can also be unlined.

Spinal Reliefs

Surgical relief (wide and gradual)

Spinous Process (narrow)

Fasteners

Screws and Barrel Nuts

Speedy Rivets

Other available options

Paraspinal reinforcements