



Anterior View



Left Lateral
View



Posterior View

Product Overview

The LSO Flex Foam® 1, anterior opening with tongue, is indicated for a number of maladies ranging from spinal stenosis to degenerative disc disease and spondylolisthesis. It's efficient in trunk support for those patients who are restricted to wheelchair use.

There are multiple opening and frame designs available.

Spinal Indications

- Herniated Disc
- Low Back Pain
- Lumbar Instability
- Degenerative Disc Disease
- Spondylolysis
- Spinal Stenosis
- Compression Fracture

Features

- Patented Flex Foam® design.
- Ease of donning and doffing.
- Lightweight.
- Removable external frame allows for easy adjustment.

Send your orders to:
orders@spinaltech.com

Send your CAD files to:
cad@spinaltech.com

**For product options see
reverse side.**

Product Options

Recommended Plastic

Copolymer (CP) 1/8" or 5/32" depending on patient size (Call for available color options)

Available Plastics

Polyethylene (PE) 3/32", 1/8", 5/32", 3/16" (Least Rigid)

Modified Polyethylene (MPE) 1/8", 5/32", 3/16" (Semi Rigid) (if Transfer pattern is requested)

Standard Flex Foam® Liner

1/4" (soft) x 1/8" (firm)

Other Liner Thickness Available

3/16" (soft) x 1/8" (firm)

1/8" (soft) x 1/8" (firm)

Spinal Reliefs

Surgical relief (wide and gradual)

Spinous Process (narrow)

Fasteners

Screws and Barrel Nuts

Speedy Rivets

Send your orders to:
orders@spinaltech.com

Send your CAD files to:
cad@spinaltech.com

3.24

Flex Foam® 1 - Frame Design Options



Low Profile Frame
Trim Lines with
More Rigid Lateral
Support



Low Profile Frame
Trim Lines with
More Rigid Lateral
Support - Side View



Low Profile Frame
Trim Lines Allowing
Some Flexibility
Along with Support



Low Profile Frame
Trim Lines Allowing
Some Flexibility
Along with Support

Send your orders to:
orders@spinaltech.com

Send your CAD files to:
cad@spinaltech.com

3.24