

# Single Opening With Hip Joint and Thigh Cuff



Anterior View



Lateral View

Posterior View

## Product Overview

The orthosis design shown is a single opening LSO with a hip joint and thigh cuff. This style allows for greater adjustability while limiting motion. Adjustment bars are positioned on orthosis and thigh section for adjustability.

## Spinal Indications

- Immobilization of Lumbosacral joint
- Pre and Post-operative hip surgery
- Pre and Post-operative spinal situations
- Spinal fusion

## Features

- Choice of three hip joints available.
- Reinforced hip and thigh support
- Adjustment bars.
- Easy positioning of hip joint

Send your orders to:  
[orders@spinaltech.com](mailto:orders@spinaltech.com)

Send your CAD files to:  
[cad@spinaltech.com](mailto:cad@spinaltech.com)

**For product options see reverse side.**

## **Product Options**

### **Recommended Plastic**

Modified Polyethylene (MPE) 1/8", 5/32", 3/16" (Semi Rigid)

---

### **Available Plastics**

Copolymer (CP) 1/8" or 5/32" depending on patient size  
(Call for available color options)

---

### **Recommended Aliplast Liner**

3/16" Aliplast

---

### **Other Liner Thickness Available**

1/8", 1/4" Aliplast  
Note: This orthosis can also be unlined.

---

### **Spinal Reliefs**

Surgical Relief (wide and gradual)  
Spinous Process (narrow)

---

### **Fasteners**

Screws and Barrel Nuts  
Speedy Rivets

---

### **Available Hip Joints**

Becker Drop Lock  
Becker Multi Axis  
Other Hip Joints Available | Call us for more information

---

### **Other available options**

Paraspinal reinforcements

Send your orders to:  
[orders@spinaltech.com](mailto:orders@spinaltech.com)

---

Send your CAD files to:  
[cad@spinaltech.com](mailto:cad@spinaltech.com)

5.24