



# The Next Generation of Comfort & Effortless Adjustment

---



**Spinal  
Technology**



**Series** of  
Spinal Orthoses

# SpinaLoc®

## Overview

SpinaLoc® system of prefabricated braces is an ideal remedy when an immediate fitting is necessary. SpinaLoc® is suitable for treatment of low back pain, degenerative disc disease and lumbar instability along with lumbar compression fractures.

## Features

The two piece construction provides structural integrity. The dynamic enclosure system ensures ease of application and adjustability.


- Removable anterior and posterior panels can be heat formed to customize fit
- Rigid molded side panels for lateral stability
- Anterior section with full abdominal available
- Conversion to a TLSO can be easily achieved by adding ATE and/or PTE

## Accessories

- ATE w / Axilla Straps
- ATE w / Shoulder Straps
- Anterior Panels
- Posterior Panels



**L0637**



Adjustable Hook and Loop Straps



Rigid Posterior Panel



Hand Strap for Ease of Donning



Adjustable Hook and Loop Closure

# SpinaLoc® X1

## Overview

The SpinaLoc® X1 features comfort, stability and a variety of options. The single piece design makes for effective compression with easy donning and doffing. Pre-molded anterior inserts come neutral or pendulous and are pre-slotted to accommodate an anterior thoracic extension. The anterior insert pocket is adjustable providing correct positioning and accommodating a wider range of sizes.

The orthosis is designed with two smooth cross-pull compression straps, with elastic tension ends as well as hand-pull tabs on all strap ends.

## Features

- Posterior insert comes with spinal relief cut out.
- Anterior insert is pre-stamped to accommodate a sternal kit.
- Two smooth cross-pull compression straps with elastic tension ends make for easy application and effective compression.
- Hand-pull tabs on all ends.



**L0637  
or L0650**


### Model X1

LSO  
9" Front and 14" Back  
with lateral panels

## Accessories

### Anterior Panel Options Neutral or Pendulous

7" Neutral & Pendulous Front



### Posterior Panel Options

3 Degrees of Lordosis 5°, 15°, 25°



# SpinaLoc® Size Chart

## SpinaLoc®

Back

Front

Female Size	Hip, Waist & Xiphoid Circumference	Anterior Height	Posterior Height
F1 (S/M)	25" - 37"	9.5"	15"
F2 (M/L)	33" - 45"	10.5"	15"
F3 (L/XL)	41" - 53"	11.5"	16"



ATE with Axilla Straps



ATE / PTE with Shoulder Straps

Male Size	Hip, Waist & Xiphoid Circumference	Anterior Height	Posterior Height
M1 (S)	26" - 40"	11"	16"
M2 (M)	32" - 46"	12"	16"
M3 (L)	40" - 54"	13"	16"
M4 (XL)	44" - 58"	13"	18"



PTE with Axilla Straps

## SpinaLoc® X1 Models

Size	HIPS
Small	32" - 38"
Medium	36" - 42"
Large	40" - 46"
X-Large	44" - 50"
2X-Large	48" - 54"
3X-Large	52" - 58"
4X-Large	56" - 62"