

Pt Name:

Order #:

Ordered By:

Cust. Acct. #:

For finish trim measurements, full-time scoliosis orthoses



Standing

Attach waist band, and tighten to clearly delineate the waist line



Do not loosen or remove waist band. This will be the waist reference for all heights.

Note: Waist Belt is positioned above the Iliac Crest and below floating ribs 12 & 11.

Note: Keep in mind that superior heights will increase if orthosis is correcting scoliotic curve. This will elongate the spine.



Sitting



Extend Measuring Tool to surface of the "Firm" Chair