

Custom Spinal TLSO

Anterior Opening Orthosis

With Tongue



Anterior View



Left Lateral
View



Posterior View

Product Overview

This TLSO utilizes a single anterior opening. One piece of plastic is measured to wrap around the patient and tighten in the front (anterior), using three straps. An additional piece of plastic (tongue) is placed behind the straps.

Spinal Indications

- Post-Operative Lumbar Surgery
- Degenerative Disc Disease
- Spondylolysis
- Spondylolisthesis
- Spinal Stenosis
- Compression Fracture
- Herniated Disc
- Low Back Pain

Features

- More total circumferential control.
- Ease of donning and doffing.

Opening Options

Anterior Opening with Tongue (Pictured), Posterior, Lateral (Right or Left)

Send your orders to:
orders@spinaltech.com

Send your CAD files to:
cad@spinaltech.com

For product options see
reverse side.

Product Options

Recommended Plastic

Polyethylene (PE) 3/32", 1/8", 5/32", 3/16" (Least Rigid)
Modified Polyethylene (MPE) 1/8", 5/32", 3/16" (Semi Rigid) (if Transfer pattern is requested)

Available Plastics

Copolymer (CP) 1/8" or 5/32" depending on patient size
(Call for available color options)

Recommended Aliplast Liner

3/16" Aliplast

Other Liner Thickness Available

1/8", 1/4" Aliplast
Note: This orthosis can also be unlined.

Spinal Reliefs

Surgical relief (wide and gradual)
Spinous Process (narrow)

Loops & Chafes

Plastic Loop and Chafe
Extended Plastic Chafe
Extended Dacron Chafe

Fasteners

Screws and Barrel Nuts
Speedy Rivets

Other available options

Shoulder straps
Paraspinal reinforcements

Send your orders to:
orders@spinaltech.com

Send your CAD files to:
cad@spinaltech.com

ST 9-20