

# Spinal Technology Prosthetics

## TRANSFEMORAL SOCKET TO MEASUREMENT FORM

FACILITY: \_\_\_\_\_  
ADDRESS: \_\_\_\_\_  
CITY: \_\_\_\_\_ STATE: \_\_\_\_\_  
ZIP: \_\_\_\_\_ PHONE: \_\_\_\_/\_\_\_\_/\_\_\_\_  
PRACTITIONER: \_\_\_\_\_  
E-MAIL: \_\_\_\_\_

PATIENT'S NAME: \_\_\_\_\_  
 LEFT       RIGHT       BILATERAL  
AGE: \_\_\_\_\_ SEX: \_\_\_\_\_ HGT: \_\_\_\_\_ WGT: \_\_\_\_\_  
COLOR: \_\_\_\_\_ SWATCH #: \_\_\_\_\_  
ACTIVITY LEVEL: \_\_\_\_\_

P.O. #: \_\_\_\_\_ ORDER DATE: \_\_\_\_\_ DATE NEEDED: \_\_\_\_\_ SHIP METHOD: \_\_\_\_\_

### FABRICATION

#### (1) SOCKET DESIGN (from Template Brochure)

- GFH (SABOLICH STYLE)
- AGGRESSIVE IC
- MODERATE IC
- MILD IC
- TRACER CAT-CAM
- SOFT NARROW M-L
- NARROW M-L #ONE
- NARROW M-L #TWO
- NARROW M-L #THREE
- NARROW M-L #FOUR
- NARROW M-L #FIVE
- SOFT QUAD
- QUAD #ONE
- QUAD #TWO
- QUAD #THREE
- QUAD #FOUR

#### (2) SOCKET MATERIAL

- Vivak
- Thermolyn Stiff
- Polypropylene

#### (3) SUSPENSION DESIGN

- Prep for Bulldog lock (mark button side)
- Prep for Bulldog Lanyard
- Prep for KISS Lanyard
- AK Valve Prep/Install: (PA0002)
- Other:

#### (4) DISTAL ATTACHMENT

- Base plate (FHLA-6)
- Three prong
- Four prong
- Other:

#### Componentry Disclaimer:

We recommend that you provide components. We can prep the socket to accept drop-in components. If we provide any items, we will bill the item at cost plus a 15% handling fee.

### ADDITIONAL INSTRUCTIONS

Submit orders to [ProstheticOrders@SpinalTech.com](mailto:ProstheticOrders@SpinalTech.com)

Spinal.Tech/Prosthetics

Follow Us




8075 National Turnpike,  
Louisville, KY 40214  
(855) 569-3754

## PATIENT MEASUREMENTS

ANTERIOR

LATERAL

SOUND SIDE



**DISTAL END CONTOUR:**

- ROUND
- FLAT
- CONICAL

- A. LENGTH TO DISTAL FEMUR
- B. LENGTH TO DISTAL END (UNCOMPRESSED)
- C. SOFT TISSUE M-L (BELOW ISCHIUM)
- D. BONY M-L (RAMUS TO LATERAL WALL)
- E. A-P
- F. GLUTEAL FOLD TO DISTAL END
- G. ISCHIUM TO FLOOR

- H. KNEE CENTER TO FLOOR
- I. PTB TO FLOOR
- J. CALF CIRCUMFERENCE
- K. HEIGHT OF CALF
- L. ANKLE CIRCUMFERENCE
- M. HEIGHT OF ANKLE

### ADDITIONAL INSTRUCTIONS

Submit orders to [ProstheticOrders@SpinalTech.com](mailto:ProstheticOrders@SpinalTech.com)

Spinal.Tech/Prosthetics

Follow Us



**8075 National Turnpike,  
Louisville, KY 40214  
(855) 569-3754**