



**Spinal
Technology**
Prosthetics

CAD FABRICATION BROCHURE ETS TEMPLATES

Prosthetic Central Fabrication

Electronic Test Sockets

Clinical & Technical Support

Spinal.Tech/Prosthetics

Follow Us    

Submit Orders to
ProstheticOrders@SpinalTech.com

**8075 National Turnpike,
Louisville, KY 40214**

(855) 569-3754



**Spinal
Technology**
Prosthetics

Spinal Technology is proud to announce our new Prosthetics division. **Glenn Hutnick and his team** of experienced central fabrication prosthetic technicians have moved to our Louisville, Kentucky location to continue manufacturing prosthetic sockets.

Together, our ABC-certified staff of Orthotists and Prosthetists collaborate with highly skilled, experienced technicians to provide the highest quality custom products and prosthetic devices, fast delivery time, and unparalleled customer support.

Using **computer-aided design (CAD)** modification and fabrication, our experienced team has been making prosthetic sockets for over 25 years. We specialize in both above-the-knee (AK) and below-the-knee (BK) diagnostic and definitive sockets. All orders are custom fabricated to your patient's cast, scans, and measurements.



**Spinal
Technology**
Prosthetics

Your Full Service Prosthetic Fabricator

CAD SOCKET TEMPLATES

The following pages contain images of various socket shapes and designs that we can produce for your patient from measurements.

Please be aware that we can also image any socket shape that you may prefer and library it for your own use. We can then provide that shape for you using the measurements you provide.

We also have the capacity to modify the mold from measurement to correct any particular problems your patient may have. Additional relief on the lateral distal end can be provided, or hugging in the lateral wall. We can even accommodate flexion and abduction problems.

All we need from you is a completed form that describes, in detail, what you would like.

Please feel free to call us to discuss these and other central fabrication services that we can support you and your patient with.

SOCKET SHAPES AVAILABLE

- | | |
|--------------------------|----------------------|
| 1) GFH MILD IC | 9) NARROW M-L THREE |
| 2) GFH MODERATE IC | 10) NARROW M-L FOUR |
| 3) GFH AGGRESSIVE IC | 11) NARROW M-L FIVE |
| 4) GFH CAT-CAM DESIGN | 12) SOFT QUAD DESIGN |
| 5) TRACER CAT-CAM DESIGN | 13) QUAD ONE |
| 6) SOFT NARROW M-L | 14) QUAD TWO |
| 7) NARROW M-L ONE | 15) QUAD THREE |
| 8) NARROW M-L TWO | 16) QUAD FOUR |



Spinal.Tech/Prosthetics

Follow Us    


**8075 National Turnpike,
Louisville, KY 40214**

(855) 569-3754


MILD ISCHIAL CONTAINMENT SOCKET TEMPLATE




A



L





M




P


MODERATE ISCHIAL CONTAINMENT SOCKET TEMPLATE




A



L





M




P


AGGRESSIVE ISCHIAL CONTAINMENT SOCKET TEMPLATE




A



L




M




P


GFH CAT-CAM DESIGN




A



L




M



P

TRACER CAT-CAM DESIGN



A



L




M




P

SOFT NARROW M-L




A



L




M



P

NARROW M-L ONE



A



L




M



P

NARROW M-L TWO



A



L




M



P

NARROW M-L THREE



A



L




M



P

NARROW M-L FOUR



A



L




M




P

NARROW M-L FIVE



A



L




M



P

SOFT QUAD DESIGN



A



L




M




P

QUAD ONE



A



L





M




P


QUAD TWO




A



L




M




P


QUAD THREE



A



L




M



P


QUAD FOUR



A



L



M



P



**Spinal
Technology**
Prosthetics



***Glenn Hutnick,
CPO, CTP***

Glenn has been a Certified Practitioner since 1979. He opened his clinical practice in 1980 and has focused his career on prosthetics in socket development and fabrication methods. Glenn was the founding father of NAAOP in 1987 and its first president.

Glenn and his team bring over 40 years of clinical, technical, research, and education experience in Prosthetics and Orthotics, including 25 years of CAD modification and fabrication experience to the team at Spinal Technology.

When Glenn is not working, he enjoys golf, fishing, traveling, and cooking with his wonderful wife Diane, two adult sons Jonathan and Gregory, and daughter-in-law Alicia.



Scan the QR code
to find all our prosthetic order forms.

Spinal.Tech/Prosthetics

Follow Us

Submit Orders to
ProstheticOrders@SpinalTech.com

8075 National Turnpike,
Louisville, KY 40214

(855) 569-3754